

Exercise 2: Favicons

All work in this exercise needs to happen on your USB flashdrive.

1. In TextWrangler, open index.html.
2. In a browser, go to this URL: https://www.sugarsync.com/pf/D242105_6103826_8384214.
3. Download, unzip the file and copy it into your images folder.
4. In Photoshop, crop the image to a square and save it in PNG format.
5. In your browser, go to the Favicon Generator <http://tools.dynamicdrive.com/favicon/>.
6. Upload your PNG image, and generate a favicon. Download the favicon and place it at the root level of your web folder.
7. At the bottom of the page in the Favicon Generator, there is one line of code: copy and paste it into the <head> section of index.html.
8. Save the file and load the page in your browser. You should see it in the tab of your current window.